

stellarex™

Drug-coated
Angioplasty Balloon

Stellarex™

0.014" OTW Drug-coated Angioplasty Balloon

**Built on a Clinically Proven
SFA Coating Technology.**

**The
Clear DCB
Choice**

*0,36 mm

Spectranetics®
Always Reaching Farther

Stellarex™ 0.014"* OTW Drug-coated Angioplasty Balloon

Product Catalog Number	Sheath Size (Fr)	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
AA14SX020040090	4	2	40	90	8	14
AA14SX020080090	4	2	80	90	8	14
AA14SX020150090	4	2	150	90	8	14
AA14SX025040090	4	2.5	40	90	8	14
AA14SX025080090	4	2.5	80	90	8	14
AA14SX025150090	4	2.5	150	90	8	14
AA14SX030040090	4	3	40	90	8	14
AA14SX030080090	4	3	80	90	8	14
AA14SX030150090	4	3	150	90	8	14
AA14SX035040090	4	3.5	40	90	8	14
AA14SX035080090	4	3.5	80	90	8	14
AA14SX035150090	4	3.5	150	90	8	14
AA14SX040040090	4	4	40	90	8	14
AA14SX040080090	4	4	80	90	8	14
AA14SX040150090	4	4	150	90	8	14
AA14SX020040150	4	2	40	150	8	14
AA14SX020080150	4	2	80	150	8	14
AA14SX020150150	4	2	150	150	8	14
AA14SX025040150	4	2.5	40	150	8	14
AA14SX025080150	4	2.5	80	150	8	14
AA14SX025150150	4	2.5	150	150	8	14
AA14SX030040150	4	3	40	150	8	14
AA14SX030080150	4	3	80	150	8	14
AA14SX030150150	4	3	150	150	8	14
AA14SX035040150	4	3.5	40	150	8	14
AA14SX035080150	4	3.5	80	150	8	14
AA14SX035150150	4	3.5	150	150	8	14
AA14SX040040150	4	4	40	150	8	14
AA14SX040080150	4	4	80	150	8	14
AA14SX040150150	4	4	150	150	8	14

*0,36 mm

Indications

The Stellarex™ OTW drug-coated angioplasty balloon is indicated for the treatment of de-novo or re-stenotic lesions in the lower extremities to establish blood flow and to maintain vessel patency.

Important Safety Information

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Corporate Headquarter

The Spectranetics Corporation
9965 Federal Drive
Colorado Springs, CO 80921
Tel: +1 719-447-2000

Spectranetics International BV

Plesmanstraat 6
3833 LA Leusden
The Netherlands
Tel: +31 33 43 47 050

Spectranetics France S.A.R.L.

18 rue de l'Arcade
75008 Paris
France
Tel: +33 1 40 98 08 73

Spectranetics Deutschland GmbH

Schweinfurter Strasse 7
97080 Würzburg
Germany
Tel: +49 931.4520080

Learn more and come visit us at www.stellarexdcb.com

©2017 Spectranetics. All rights reserved. Approved for external Distribution. D035969-02 082017

Material not intended for distribution in France and the U.S.

Stellarex™ is a Class III medical device, CE0086.